

# Invacare® | Solara® 3G



## Tilt-In-Space Wheelchair

The *Invacare® Solara®* Wheelchair is known for its durability of frame construction and gear tilt mechanism design. The new generation *Invacare® Solara® 3G* Wheelchair retains this durability, and once again revolutionizes tilt-in-space wheelchair design with incredible simplicity, portability and ease of use and adjustment.

**Code:** SOLARA3G



## Standard features

- 136kg standard, 180kg weight capacity (heavy-duty package)
- Center-of-gravity remains constant throughout tilt range
- Positive lock mechanism eliminates “sponginess” common to locking hydraulic cylinders
- Tilt angle indicator
- Standard flared front provides more space for user’s feet and reduces interference when tilting in the anterior direction
- One axle plate accommodates all wheel sizes
- -5° to 50° of total tilt available
- Adjustable-depth, width and adjustable angle back canes



## Options



The **DeGage Dynamic Rocker Back** accommodates more aggressive rockers, like children with uncontrolled extension of the upper and lower extremities. The Rocker Back will allow an individual to move in the direction they are pushing and as a result will help to reduce pain or discomfort. For more information please see [www.degage.us](http://www.degage.us).

## Technical Data

### Limited Warranty:

- 5 years frame
- 1 year other components



30 - 60 cm



30 - 55 cm



5° - +50°



43 - 60 cm



Seat width: + 21 cm



136kg standard  
180kg wHD package\*\*



From 16 kg

\*\* Heavy duty package provides reinforced frame.  
All dimensions are also available without heavy duty package.

## Frame Colour

Please note colors may vary slightly from those displayed here.



Electric Blue



Grasshopper Green



Sunny Yellow



Electric Red



Bubblegum Pink



Grape Madness



Black Prism



Wet Black



Midnight Blue



Silver Metallic



White Snow



Tangerine Red



Sunkissed Orange