

rea | *Azalea / Azalea Assist*



Offers excellent seating comfort and many adjustment possibilities

The *Rea Azalea* is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes *Azalea* the ideal choice. *Azalea* is a self propelling wheelchair with large rear wheels, and is complimented with an alternative transit version that keeps the same frame but with small rear wheels called *Rea Azalea Assist*. The new backrest, *Flex 3*, is multi adjustable in height, width and shape for different body forms.



Codes: AZALEA-SP-XX (self propel)
AZALEA-TR-XX (transit)
(refer to customer service for the code size you require)



Features and options

Tilt & recline makes life comfortable

The *Rea Azalea* and *Rea Azalea Assist* have 25° seat tilt and a 30° backrest recline. Both tilt and recline can be easily adjusted manually or electrically by the assistant and also by the user for the electrical version.

The backrest has a biomechanical positioning hinge, which follows the user's natural back movement when adjusting the backrest angle. Combined with the seat-tilt function, the DSS (Dual Stability System) helps to prevent Decubitus thanks to its nonfriction re-body positioning.



Flex 3, backrest plate
Apart from being height adjustable it is now width adjustable as well.



Flex 3, tension adjustable backrest
The multi adjustable backrest can be adjusted in height and in width for different body shapes.



Wide range of options
Individual width adjustments in the armrests and legrests and accessories for individual options.



Power tilt
The hand control for tilt and recline can be used by the attendant or user.

Technical data

Rea Azalea & Azalea Assist		390-450/440-500/ 490-550 mm (390-490/440-540/ 490-590 mm)		430 - 500 mm		-1° to 23° (+1° to 25°)		400/450 mm (450 mm)		0°-30°		Backrest plate: 560-790 mm Tension adjustable backrest: 540-715 mm		330 - 500 mm
		240 - 360 mm		SW + 220 mm (SW + 230 mm)		950-1020 mm		900-1250 mm		Max 135 kg		34 kg (SW 440 mm)		Transport weight 16.5 kg

The Rea Azalea and Azalea Assist complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.

Frame colour

Please note colors may vary slightly from those displayed here.



Grey

Upholstery



Black plush TR35