

Invacare® **TDX SP**



A quick and quiet superior ride

The *Invacare*® **TDX**® SP wheelchair is the first model in the next generation of TDX (Total Driving eXperience) power wheelchairs. The TDX SP wheelchair takes all the best features of the TDX series and improves upon them with enhanced performance, superior ride quality, quieter chassis and an elegantly simple design.

By optimising the Torque speed curve, the motor/electronics combination provides a 9.3km/h top speed with ample torque. The chair gets consumers where they need to go quickly, while also providing the power required to handle those tough obstacles.

Seating Options

The *Invacare® TDX® SP* wheelchair offers a wide variety of seating options that make it an extremely versatile power chair. Various power seating options, refer to the **Invacare Powered Seating and Positioning Systems Options Portfolio**.

- Various seat-to-floor height is as low as 42.5 - 49cm (16.75" - 19.25")
- Up to 50° of tilt
- Up to 170° of recline
- Up to 30cm (12") of elevation
- All powered seating functions can be operated either through or separate from the driver control.



Adjustable Rehab Seat



Tilt



Recline



Tilt/Recline

Also available - Tilt/Recline/Elevate, Elevate 12" and 6".

Technical Data

Specifications are approximate and may vary with options selected.



111.5 cm L**
68 cm W**
90 cm H**




28 - 61 x
25.5-61cm




42.5 - 49 cm



35.5 cm



7.5 cm



52 cm



9°



Up to 181 kg



136 kg
181 kg HD



4 Pole



9 kph
8 kph HD



Group 24,
22NF



Up to 32 kms




5 yrs Base Frame
3 yrs Seat Frame
1 yr Electronics/
Motor/Gearbox

**Rehab seat and 70° swingaway footrest

Core Technologies



Invacare® SureStep® Suspension with Stability Lock provides a smooth and stable transition



Invacare® G-Trac™ Technology provides precise tracking for imprecise environments



Invacare® MK6i™ Electronics are the benchmark in innovative technology



Motion Concepts power positioning systems can be interfaced with the Invacare TDX SP power base.