



Yes, you can.®

Price List and Order Form

Effective Feb, 2009

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group _____ Date _____

Therapist/Email _____ Therapist Phone _____

Client _____ Therapist Fax _____

Ship to Name & Address _____

POWERED SEATING BASE

- 3GTQ Torque SP with Formula CG Tilt Only 10920.00
- 3GTQ Torque SP with Formula CG Elevate Only 10920.00
- 3GTQ Torque SP with Formula CG Tilt/Recline..... 13030.00

- W17 17" wide.....No Charge
- W18 18" wide.....No Charge
- W19 19" wide.....No Charge
- W20 20" wide.....No Charge
- W21 21" wide.....No Charge
- W22 22" wide.....No Charge
- W23 23" wide.....No Charge
- W24 24" wide.....No Charge

Base Package

- 3GTQ-CG Torque SP for Formula CG Tilt in Space with MK6i
- 136 kg Weight Capacity
- Performance 4 Pole MotorsStandard
- 8" X 2.25" Semi Pneumatic CastersStandard
- Flip Back Adjustable Height ArmStandard
- MK690ACC Controller w/Actuator ControlStandard
- 70NHD 70° Swingaway FootrestsStandard
- 8 Amp Off Board Battery ChargerStandard
- 1515 Push Button Style Seat Positioning StrapStandard
- Crash Tested Transport BracketsStandard

ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Deep (Adjustable).....No Charge
- 1619D 16" to 19" Deep (Adjustable)..... Standard
- 1922D 19" to 22" Deep (Adjustable).....No Charge

ADJUSTABLE SEAT DEPTH SETTING

- D12 Seat Depth 12" (Short Base).....No Charge
- D13 Seat Depth 13" (Short Base).....No Charge
- D14 Seat Depth 14" (Short Base).....No Charge
- D15 Seat Depth 15" (Short Base).....No Charge
- D16 Seat Depth 16" (Short Base).....No Charge
- D17 Seat Depth 17" (Short Base).....No Charge
- D18 Seat Depth 18" (Short Base).....No Charge
- D19 Seat Depth 19".....No Charge
- D20 Seat Depth 20".....No Charge
- D21 Seat Depth 21".....No Charge
- D22 Seat Depth 22".....No Charge

HEAVY DUTY OPTION 181 KG CAPACITY

- 3GRX-CG Ranger X – Tilt only – Heavy Duty Motors..... 1,650.00

SYSTEM TYPE

- CGT CG Tilt Only SystemStandard
- CNVACT Conventional Actuator.....Standard
- SMTACT Smart Actuator..... 715.00

SEAT-TO-FLOOR HEIGHT

- LOW 17.75" at 0°.....No Charge
- MED 18.75" at 0°.....Standard
- TALL 19.75" at 0°.....No Charge

SEAT TILT⁵

- OSA 0° - 55°..... Standard
- SSA 5° - 60°.....No Charge

SEAT SELECTION

- BPO Back Posts Only – Round Canes Only.....No Charge
- CURVABS 435.00

BACK HEIGHT

- TSBH16 Back Height 16".....No Charge
- TSBH17 Back Height 17".....No Charge
- TSBH18 Back Height 18".....No Charge
- TSBH19 Back Height 19".....No Charge
- TSBH20 Back Height 20"..... Standard
- TSBH21 Back Height 21".....No Charge
- TSBH22 Back Height 22".....No Charge
- TSBH23 Back Height 23".....No Charge
- TSBH24 Back Height 24".....No Charge
- TSBH25 Back Height 25".....No Charge
- TSBH26 Back Height 26".....No Charge

SEAT WIDTH RANGE

- 1216W 12" to 16" Width (Adjustable).....No Charge
- 1620W 16" to 20" Width (Adjustable).....Standard
- 2024W 20" to 24" Width (Adjustable).....No Charge

ADJUSTABLE SEAT WIDTH SETTING

- W12 12" wide.....No Charge
- W13 13" wide.....No Charge
- W14 14" wide.....No Charge
- W15 15" wide.....No Charge
- W16 16" wide.....No Charge

FOOTNOTES

- ¹ Not available with Transport Tie Down option.
- ² Push handles included with BPO.
- ³ Seat Widths 12" – 15" and Seat Depths 12" – 15" are only available with Back Posts Only or Conventional Seating. For users under 150 lbs. These sizes are not available with Contoura Seating.

FOOTNOTES

- ⁴ Available only with 181 kg package.
- ⁵ Indicates static tilt of seat, impacts seat to floor height. Every 5 degrees raises the seat to floor height.
- ⁶ For Curved ABS only.
- ⁷ Not available with seat widths 12" – 15".
- ⁸ For Conventional back only.

BASE

CASTER OPTIONS

- 1259 6" x 2" Semi Pneumatic Casters No Charge
- 1228-3 8" x 2.25" Semi Pneumatic Casters Lt. Grey Standard
- 1251-3 8" x 2" Pneumatic Casters No Charge
- 1220-3 9" x 2.75" Pneumatic Casters No Charge ⁹
- 1465-3 9" x 2.75" Foam Filled Casters No Charge ⁹

CASTER SUSPENSION ¹⁰

- 9213-8 Shock Fork for 8" Casters 149.00

WHEEL OPTIONS

- 1430-3 14" x 3" Pneumatic Tires Standard
- 1431-3 14" x 3" Tires w/Foam Filled Inserts 275.00 ¹¹
- 1440 14" x 4" Pneumatic Tires 220.00 ¹¹
- 1441 14" x 4" Tires w/Foam Filled Inserts 352.00 ¹¹

FRAME FINISH

- 61P Electric Red Standard
- 62P Electric Blue Standard

HEADREST TYPE

STEALTH COMFORT HEADRESTS

- CP175 14" Multi Axis Headrest 250.00
- CP175HW Removable & Adjustable Hardware 250.00
- CP275 10" Multi Axis Headrest 250.00
- CP275HW Removable & Adjustable Hardware 250.00

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj. Hgt No Charge
- 29F Desk Length Adj. Hgt No Charge

ELECTRONICS

EXPANDABLE CONTROLLER SYSTEMS ¹²

- MPJM6 MPJ+ Multiple Drive Joystick Standard
- ASLRDYM6MK6 Display Only (No Driver Control) 825.00 ¹³

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ 435.00 ¹⁴
- LEFT Left Hand Mounted Joystick No Charge

POWERED SEATING ACCESSORIES

- SANODE Single Act. Oper. through Driver Ctrl Standard ¹⁵
- S4WSB Multiple Actuator Interface Box - must be Selected if more than one actuator function i.e. Tilt and Elevate Seat 715.00
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge

FOOTNOTES

⁹ Not available with short base.

¹⁰ Not available with 9" casters.

¹¹ 14" x 4" wheels increase the overall width of the base by 2.5". Does not include 3G design wheel or fender. **Not available with (TRUTK) TrueTrack Package.**

¹² Includes Monoport "Y" Cable to provide additional switch input port.

¹³ Display Only, No Driver control.

¹⁴ QLA and QLAM6 include ARM250 standard.

¹⁵ Required for Tilt Operation through Driver Control – Single Actuator Systems Only. Not compatible with SPJ+ joysticks. Does not come standard with egg switch.

FRONT RIGGINGS

FRONT RIGGING

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests Standard ¹⁶
- 70STAPER70° Tapered Footrests 127.00
- AT5544 Manual Elevating Legrests 352.00
- Stump Support Right 215.00
- Stump Support Left 215.00

POWER ELEVATING LEGRESTS-ARTICULATING

- ELRPW Power ELR w/Articulation 20" Ext 2195.00 ¹⁷

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip-Up Composite Footplates Standard
- AT5543 Adjustable Angle Flip-Up Footplates 205.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf 39.00
- 1600BK Heel Loops w/Ankle Straps 72.00
- IMPACTG Impact Guards 72.00
- ALPT4 Longer Pivot And Slide Tube 4" 99.00

MISCELLANEOUS

MISCELLANEOUS ACCESSORIES

- 2ATP Perspex Tray – Swing away 545.00

FOOTNOTES

¹⁶ Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.

¹⁷ Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.