

# Invacare® 3G Torque Tilt in Space - Category 4 RWD TNS



Yes, you can.®

## Price List and Order Form

Effective Feb, 2009

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group \_\_\_\_\_ Date \_\_\_\_\_

Therapist/Email \_\_\_\_\_ Therapist Phone \_\_\_\_\_

Client \_\_\_\_\_ Therapist Fax \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

### POWERED SEATING BA SE

- 3GTQ-CGTorque SP for Formula CG with Tilt MK6i ..... 10920.00

Base Package	
>	3GTQ-CG Torque SP for Formula CG Tilt in Space with MK6i
>	136 kg Weight Capacity
>	Performance 4 Pole Motors ..... Standard
>	8" X 2.25" Semi Pneumatic Casters ..... Standard
>	Flip Back Adjustable Height Arm ..... Standard
>	MK690ACC Controller w/Actuator Control ..... Standard
>	70NHD 70° Swingaway Footrests ..... Standard
>	8 Amp Off Board Battery Charger ..... Standard
>	1515 Push Button Style Seat Positioning Strap ..... Standard
>	Crash Tested Transport Brackets ..... Standard

### HEAVY DUTY OPTION 181 KG CAPACITY

- 3GRX-CG Ranger X – Tilt only – Heavy Duty Motors..... 1,650.00

### SEAT-TO-FLOOR HEIGHT

- LOW 17.75" at 0° ..... No Charge <sup>1</sup>
- MED 18.75" at 0° ..... Standard
- TALL 19.75" at 0° ..... No Charge

### SEAT SELECTION

- BPO Back Posts Only – Round Canes Only ..... No Charge <sup>2</sup>
- CURV Curved Back ..... 435.00

### SEAT WIDTH RANGE

- 1216W 12" to 16" Width (Adjustable) ..... No Charge
- 1620W 16" to 20" Width (Adjustable) ..... Standard
- 2024W 20" to 24" Width (Adjustable) ..... No Charge

### ADJUSTABLE SEAT WIDTH SETTING

- W12 12" wide ..... No Charge
- W13 13" wide ..... No Charge
- W14 14" wide ..... No Charge
- W15 15" wide ..... No Charge
- W16 16" wide ..... No Charge
- W17 17" wide ..... No Charge
- W18 18" wide ..... No Charge
- W19 19" wide ..... No Charge
- W20 20" wide ..... No Charge
- W21 21" wide ..... No Charge
- W22 22" wide ..... No Charge
- W23 23" wide ..... No Charge
- W24 24" wide ..... No Charge

### ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Deep (Adjustable)..... No Charge
- 1619D 16" to 19" Deep (Adjustable)..... Standard
- 1922D 19" to 22" Deep (Adjustable)..... No Charge

### ADJUSTABLE SEAT DEPTH SE TTING

- D12 Seat Depth 12" (Short Base) ..... No Charge
- D13 Seat Depth 13" (Short Base) ..... No Charge
- D14 Seat Depth 14" (Short Base) ..... No Charge
- D15 Seat Depth 15" (Short Base) ..... No Charge
- D16 Seat Depth 16" (Short Base) ..... No Charge
- D17 Seat Depth 17" (Short Base) ..... No Charge
- D18 Seat Depth 18" (Short Base) ..... No Charge
- D19 Seat Depth 19" ..... No Charge
- D20 Seat Depth 20" ..... No Charge
- D21 Seat Depth 21" ..... No Charge
- D22 Seat Depth 22" ..... No Charge

### SEAT TILT

- OSA 0° - 55° ..... Standard
- SSA 5° - 60° ..... No Charge

### BACK HEIGHT

- TSBH16 Back Height 16" ..... No Charge
- TSBH17 Back Height 17" ..... No Charge
- TSBH18 Back Height 18" ..... No Charge
- TSBH19 Back Height 19" ..... No Charge
- TSBH20 Back Height 20" ..... No Charge
- TSBH21 Back Height 21" ..... No Charge
- TSBH22 Back Height 22" ..... No Charge
- TSBH23 Back Height 23" ..... No Charge
- TSBH24 Back Height 24" ..... No Charge
- TSBH25 Back Height 25" ..... No Charge
- TSBH26 Back Height 26" ..... No Charge

### CASTER OPTIONS

- 1228-3 8" x 2.25" Semi Pneumatic Casters Lt. Grey ..... Standard
- Other caster sizes available on request.

### WHEEL OPTIONS

- 1430-3 14" x 3" Pneumatic Tires ..... Standard

### FRAME FINISH

- 61P Electric Red ..... Standard
- 62P Electric Blue ..... Standard

### FOOTNOTES

<sup>1</sup> Not available with Transport Tie Down option.

<sup>2</sup> Push handles included with BPO.

## HEADREST TYPE

### STEALTH HEADRESTS

- CP175 14" Multi Axis Headrest ..... 250.00
- CP175HW Removable & Adjustable Hardware ..... 250.00
- CP275 10" Multi Axis Headrest ..... 250.00
- CP275HW Removable & Adjustable Hardware ..... 250.00

### ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj. Hgt ..... No Charge
- 29F Desk Length Adj. Hgt ..... No Charge

## ELECTRONICS

### EXPANDABLE CONTROLLER SYSTEMS

- MPJM6 MPJ+ Multiple Drive Joystick ..... Standard

### JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ ..... 435.00
- LEFT Left Hand Mounted Joystick ..... No Charge

## FRONT RIGGINGS

### FRONT RIGGING

#### CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests ..... Standard
- AT5544 Manual Elevating Legrests ..... 350.00
- Stump Support Right ..... 215.00
- Stump Support Left ..... 215.00

#### FOOTPLATE OPTIONS (PAIR)

- 1651 Flip-Up Composite Footplates ..... Standard
- AT5543 Adjustable Angle Flip-Up Footplates ..... 205.00

#### POWER ELEVATING LEGRESTS-ARTICULATING

- ELRPW Power ELR w/Articulation 20" Ext ..... 2195.00 <sup>3</sup>

#### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... 39.00
- 1600BK Heel Loops w/Ankle Straps ..... 72.00
- IMPACTG Impact Guards ..... 72.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... 99.00

## MISCELLANEOUS

### MISCELLANEOUS ACCESSORIES

- 2ATP Perspex Tray – Swing away ..... 545.00

## FOOTNOTES

<sup>3</sup> Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.