



Yes, you can.®

Price List and Order Form

Effective Feb, 2009

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group _____ Date _____

Therapist/Email _____ Therapist Phone _____

Client _____ Therapist Fax _____

Ship to Name & Address _____

POWERED SEATING BASE

- TDXSP TDXSP with Formula CG Tilt w/ MK6i™ 14135.00
- CGTE TDXSP CG Elevate to combine with Tilt..... 16330.00
- CGTR TDXSP CG Recline to combine with Tilt 16330.00
- CGTRE TDXSP CG Elevate/Recline with Tilt 18525.00

SEAT SELECTION

BACKS

- BPO Back Posts Only.....No Charge
- CURV Curved Back 395.00

SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range (Adjustable)No Charge
- 1620W 16" to 20" Seat Width Range (Adjustable) Standard
- 2024W 20" to 24" Seat Width Range (Adjustable)No Charge

SEAT WIDTH SETTING

- W12 12" Seat WidthNo Charge
- W13 13" Seat WidthNo Charge
- W14 14" Seat WidthNo Charge
- W15 15" Seat WidthNo Charge
- W16 16" Seat WidthNo Charge
- W17 17" Seat WidthNo Charge
- W18 18" Seat WidthNo Charge
- W19 19" Seat WidthNo Charge
- W20 20" Seat WidthNo Charge
- W21 21" Seat WidthNo Charge
- W22 22" Seat WidthNo Charge
- W23 23" Seat WidthNo Charge
- W24 24" Seat WidthNo Charge

SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range (Adjustable).....No Charge
- 1619D 16" to 19" Seat Depth Range (Adjustable)..... Standard
- 1922D 19" to 22" Seat Depth Range (Adjustable).....No Charge

SEAT DEPTH SETTING

- D1212" Seat DepthNo Charge
- D1313" Seat DepthNo Charge
- D1414" Seat DepthNo Charge
- D1515" Seat DepthNo Charge
- D1616" Seat DepthNo Charge
- D1717" Seat DepthNo Charge
- D1818" Seat DepthNo Charge
- D1919" Seat DepthNo Charge
- D2020" Seat DepthNo Charge
- D2121" Seat DepthNo Charge
- D2222" Seat DepthNo Charge

Base Price Includes

- TDXSP for Formula CG Tilt with MK6i
- 136 Kg Weight Capacity (181kg Heavy Duty Option)
- 14" X 3" Pneumatic Tires Standard
- 4 Pole Motors Standard
- Flip Back Height Adjustable Arms..... Standard
- MK690ACC Controller w/ Actuator Control..... Standard
- 70NHD 70° Swingaway Footrests Standard
- MPJ+ Joystick Standard
- 1515 Push Button Style Seat Positioning Strap Standard
- Wheelchair Tie down Brackets Standard

SYSTEM TYPE

- CNVACTS Conventional Actuator Standard
- SMTACTS Smart Actuator..... 715.00

BASE

USER WEIGHT LIMIT

- U300 User Weight 136 kg Standard ¹
- U400 User Weight <181 kg (Not for Elevate/Recline) 1095.00 ²

WHEEL OPTIONS (PAIR)

- 1430-3 14" x 3" Pneumatic Tires..... Standard
- 1431-3 14" x 3" with Foam Filled Inserts 275.00

FRAME FINISH

CONTEMPORARY

- 61P Electric Red Standard
- 62P Electric Blue Standard

SEAT

SEAT TILT

- OSA 0° - 55° Seat Tilt..... Standard

FOOTNOTES

¹ Package includes 4 Pole Motors and 5.8 mph top speed with <136 kg user weight limit. Package includes Quiet Stability Lock.
² Package includes Heavy-Duty 4-pole motors and 5 mph top speed with <181 kg user weight limit.

FOOTNOTES

³ Available only with 181 kg. package.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP

Accessories and Options

BACK HEIGHT

- TSBH20 Back Height 20" Standard
- TSBH21 Back Height 21" No Charge
- TSBH22 Back Height 22" No Charge
- TSBH23 Back Height 23" No Charge
- TSBH24 Back Height 24" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH26 Back Height 26" No Charge

ADJUSTABLE ARM TYPE TILT/ELEVATE (PAIR)

- 19F Full Length Adj. Hgt. Right No Charge
- 29F Desk Length Adj. Hgt. Left Standard

ADJUSTABLE ARM TYPE WITH PWR RECLINE (PAIR)

- RA19F Reclining Full Length Adj. Hgt. Right No Charge
- RA29F Reclining Desk Length Adj. Hgt. Left Standard

STEALTH COMFORT HEADRESTS

- CP175 14" Multi Axis Headrest 250.00
- CP175HW Removable & Adj Hardware 250.00
- CP275 10" Multi Axis Headrest 250.00
- CP275HW Removable & Adj Hardware 250.00

ELECTRONICS

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6 MPJ™+ Multiple Drive Joystick Standard
- ASLRDYM6 MK6 Display Only (No Driver Control) 825.00 ⁴

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ 435.00
- LEFT Left Handed Mounted Joystick No Charge

POWERED SEATING ACCESSORIES

- SANODE Single Act. Oper. through Driver Ctrl Standard ⁵
- S4WSB Multiple Actuator Interface Box - must be Selected if more than one actuator function i.e. Tilt and Elevate Seat 715.00
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING ⁶

CENTER MOUNT ⁷

- CTMT-F Center Mount Adjustable Knee Angle 242.00 ⁸
- CTMT-FJ Jr. Center Mount Adjustable Knee to Heel 242.00 ⁹

CENTER PIVOT STYLE (PAIR) ¹⁰

- 70NHD 70° Swingaway Footrests Standard ¹¹
- 70STAPER 70° Tapered Footrests 127.00 ¹²
- AT5544 Manual Elevating Legrests 352.00 ¹³

FOOTNOTES

- ⁴ Display Only, No Driver control.
- ⁵ Required for Tilt Operation through Driver Control – Single Actuator Systems Only. Not compatible with SPJ+ joysticks. Does not come standard with egg switch.
- ⁶ All front riggings include a standard footplate .
- ⁷ Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- ⁸ Extension range 12.5"-15".
- ⁹ Extension range 9.25"-12".
- ¹⁰ Provides up to 1 ½" of depth adjustment in ½" increments.
- ¹¹ With seat widths of 21" to 24" and depth of 16", must pick footplate 1350. The AT5543 is not available..
- ¹² Standard with impact guards. Not available with widths 19" to 24" and 16" deep seat.

POWER ELEVATING LEGRESTS-ARTICULATING – INCLUDES ELECTRONICS

- MPCTMT Power Center Mount Foot Platform 2195.00 ¹⁴
 - IFP Individual Footplates 275.00 ¹⁵
- MPCTMT-L Pwr Center Mnt Foot Platform Long 2195.00 ¹⁴
- ELRPW Power ELR w/Articulation 20" Ext 2195.00 ¹⁶
- ELRPW22 Power ELR w/Articulation 22" Ext 2195.00 ¹⁷

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up Footplates Standard ¹⁸
- AT5543 Adj Angle Flip-Up Footplates 205.00 ¹⁹
- CTMTFB Footboard for Center Mount 165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap 39.00
- 1600BK Heel Loops w/Ankle Straps) 72.00
- IMPACTG Impact Guards 72.00
- ALPT4 Longer Pivot And Slide Tube 4" 99.00

MISCELLANEOUS ACCESSORIES

- 2ATP Perspex Tray – Swing away 545.00

Motion Power Center Mount Knee to Heel			
STF Pick	STF at 0 deg	KTH w/IFP on MPCTMT	KTH with Footboard
Low	17.5"	10.5" - 14.5"	12" - 13.5" MPCTMT only
Med	18.5"	11.5" - 15.5"	12" - 14.5" MPCTMT only
High	19.5"	11.5" - 16.5"	12" - 15.5" MPCTMT or MPCTMT-L

FOOTNOTES

- ¹³ Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.
- ¹⁴ Includes a standard, single footboard.
- ¹⁵ Individual footplate option for Power Center Mount. This option is only available on MPCTMT.
- ¹⁶ Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.
- ¹⁷ Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 15"-22". Sold as pair.
- ¹⁸ Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- ¹⁹ One piece footboard for Center Mount (CTMT-F) Riggings.