



Yes, you can.®

Price List and Order Form

Effective Feb, 2009

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group _____ Date _____

Therapist/Email _____ Therapist Phone _____

Client _____ Therapist Fax _____

Ship to Name & Address _____

CHAIR TYPE

- TDXSPREE TDX Spree w/5" pwr elevate w/MK6i™ 8025.00 ¹

Base Price Includes

- TDX Spree with 5" pwr elevate seat MK6i
- 75 Kg Weight Capacity
- 10.5"X 3.5" Flat Free TiresStandard
- 2 Pole Motors.....Standard
- Desk Length Height Adjustable ArmrestsStandard
- Push Button Pelvic Belt.....Standard
- MK660ACC 60 Amp ControllerStandard
- SPJ+ Single Drive JoystickStandard
- Battery Charger & 31amp BatteriesStandard
- 90NHD 90° Swingaway FootrestsStandard

BASE

SEAT-TO-FLOOR HEIGHT

- LOW 14.5"Standard
- MED 16.5"No Charge
- TALL 18.5"No Charge

FRAME FINISH

- 61P Electric Red.....Standard
- 62P Electric Blue.....Standard
- 24P Wet Black.....No Charge
- 60P Silver MetallicNo Charge
- 30P Sunny YellowNo Charge
- 119P Bubblegum PinkNo Charge
- 120P Grape Madness.....No Charge
- 62P Electric Blue.....No Charge
- 127P Grasshopper Green.....No Charge

SEAT

ADJUSTABLE SEAT WIDTH RANGE 12" – 16"

- W12 12" Seat Width.....No Charge
- W13 13" Seat Width.....No Charge
- W14 14" Seat Width.....No Charge
- W15 15" Seat Width.....No Charge
- W16 16" Seat Width.....No Charge

FOOTNOTES

¹ TDX Spree comes standard with: 75 kg user weight capacity, 8 kph top speed, 2-pole motor package & MK6i SPJ™ + Joystick, MK660ACC Controller, Quiet Stability Lock.

ADJUSTABLE SEAT DEPTH RANGE 12" – 18"

- D12 12" Seat DepthNo Charge
- D13 13" Seat DepthNo Charge
- D14 14" Seat DepthNo Charge
- D15 15" Seat DepthNo Charge
- D16 16" Seat DepthNo Charge
- D17 17" Seat DepthNo Charge
- D18 18" Seat DepthNo Charge

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj Hgt ArmNo Charge
- 29F Desk Length Adj Hgt ArmStandard

SEAT BACK TYPE

- BPO Back Posts OnlyStandard
- CURV Curved Back435.00

SEAT BACK ANGLE DEGREE

- 90BA 90°Standard
- 95BA 95°No Charge

BACK HEIGHT

- BH14 Back Height 14"No Charge
- BH15 Back Height 15"No Charge
- BH16 Back Height 16"No Charge
- BH17 Back Height 17"No Charge
- BH18 Back Height 18"Standard
- BH19 Back Height 19"No Charge
- BH20 Back Height 20"No Charge

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Pos. Strap.....Standard

ELECTRONICS

CONTROLLER SYSTEMS

- SPJAPM SPJ+ w/Joystick Actuator ControlStandard ²
- MPJM6U MPJ™+ Multiple Drive Joystick 765.00
- ASLRDYM6 MK6 Display Only (No Driver Control)825.00 ³

FOOTNOTES

² Allows actuator operation through the Driver Control. Includes port for the Attendant Switch.
³ Display Only, No Driver control. Allows choice from order form.

NOTE: All specifications and dimensions are approximate.

**Invacare® TDX® Spree Rehab Seat Power Wheelchair
Accessories and Options**

FRONT RIGGING OPTIONS

FRONT RIGGING

CENTER MOUNT

- CTMT-F Ctr Mount Adj Knee Angle (12.25"-15").....242.00 ⁴
- CTMT-FJ Ctr Mount Adj (9.25" - 12").....242.00 ⁵

CENTER PIVOT STYLE (PAIR)

- 90SWING 90° Swingaway FootrestsStandard

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up FootplatesStandard
- AT5543 Adj Angle Flip-Up Footplates205.00
- CTMTFB Footboard for Center Mount.....165.00 ⁶

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap39.00
- 1600BK Heel Loops w/Ankle Straps72.00
- IMPACTG Impact Guards.....72.00
- ALPT3 3" Extension for 90° Footrests99.00

MISCELLANEOUS ACCESSORIES

- 2ATP Perspex Tray – Swing away545.00

FOOTNOTES

- ⁴ Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Available on low seat-to-floor only.
- ⁵ Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Available on low seat-to-floor only.
- ⁶ One piece footboard for Center Mount (CTMT-F) Riggings only.

NOTE: All specifications and dimensions are approximate.