



**Yes, you can.®**

**Price List and Order Form**

Effective Feb, 2009

For ease of ordering, contact Customer Service Toll Free at: 0508-468-222 Fax: 0508-807-788

www.invacare.co.nz

Group \_\_\_\_\_ Date \_\_\_\_\_

Therapist/Email \_\_\_\_\_ Therapist Phone \_\_\_\_\_

Client \_\_\_\_\_ Therapist Fax \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

**CHAIR TYPE**

- TDXSPREE-CG TDX Spree for Formula CG w/MK6i™ ..... 9895.00 <sup>1</sup>

**Base Price Includes**

- TDX Spree for Formula CG with 5" pwr elevate seat and pwr tilt
- 75 Kg Weight Capacity
- 10.5"X 3.5" Flat Free Tires .....Standard
- 2 Pole Motors.....Standard
- Flip Back Fixed Height Arm .....Standard
- MK660ACC 60 Amp Controller .....Standard
- Battery Charger & 31amp Batteries .....Standard
- SPJPLUS SPJ+ Single Drive Joystick .....Standard
- 90NHD 90° Swingaway Footrests .....Standard
- 1515 Push Button Style Seat Positioning Strap .....Standard
- Crash Tested Transport Brackets .....Standard

**BASE**

**SEAT-TO-FLOOR HEIGHT**

- LOW 16.75" ..... Standard
- MED 17.75" ..... No Charge
- TALL 18.75" ..... No Charge

**FRAME FINISH**

- 61P Electric Red ..... Standard
- 62P Electric Blue ..... Standard
- 24P Wet Black ..... No Charge
- 60P Silver Metallic ..... No Charge
- 30P Sunny Yellow ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge
- 127P Grasshopper Green ..... No Charge

**SEAT**

**SEAT SELECTION**

- BPO Back Post Only ..... Standard <sup>2</sup>

**ADJUSTABLE SEAT WIDTH RANGE 12" – 16"**

- W12 12" Seat Width.....No Charge
- W13 13" Seat Width.....No Charge
- W14 14" Seat Width.....No Charge
- W15 15" Seat Width.....No Charge
- W16 16" Seat Width.....No Charge

**FOOTNOTES**

<sup>1</sup> TDX Spree comes standard with: 75 kg user weight capacity, 8 kph top speed, 2-pole motor package & MK6i SPJ™ + Joystick, MK660ACC Controller, Quiet Stability Lock.

<sup>2</sup> BPO standard with push handles.

**ADJUSTABLE SEAT DEPTH RANGE 12 " – 18"**

- D12 12" Seat Depth .....No Charge
- D13 13" Seat Depth .....No Charge
- D14 14" Seat Depth .....No Charge
- D15 15" Seat Depth .....No Charge
- D16 16" Seat Depth .....No Charge
- D17 17" Seat Depth .....No Charge
- D18 18" Seat Depth .....No Charge

**ADJUSTABLE ARM TYPE (PAIR)**

- 19F Full Length Adj Hgt Conv .....No Charge
- 29F Desk Length Adj Hgt Conv .....Standard

**BACK HEIGHT**

- TSBH16 Back Height 16".....No Charge
- TSBH17 Back Height 17" .....No Charge
- TSBH18 Back Height 18".....Standard
- TSBH19 Back Height 19" .....No Charge
- TSBH20 Back Height 20".....No Charge

**STEALTH COMFORT HEADRESTS**

- CP175 14" Multi Axis Headrest ..... 250.00
- CP175HW Removable & Adj Hardware ..... 250.00
- CP275 10" Multi Axis Headrest..... 250.00
- CP275HW Removable & Adj Hardware ..... 250.00

**ELECTRONICS**

**CONTROLLER SYSTEMS**

- SPJAPM SPJ+ w/Joystick Actuator Control ..... Standard <sup>3</sup>
- MPJM6U MPJ™+ Multiple Drive Joystick Incl. 16" cable. 765.00
- ASLRDYM6 MK6 Display Only (**No Driver Control**) ..... 825.00

**JOYSTICK MOUNTING HARDWARE**

- QLA Quad Link Retract. Joystick Mnt SPJ+ ..... 435.00
- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ ..... 435.00
- LEFT Left Handed Mounted Joystick .....No Charge

**FRONT RIGGING OPTIONS**

**FRONT RIGGING**

**CENTER MOUNT**

- CTMT-F Ctr Mount Adj Knee Angle (12.25"-15")..... 242.00 <sup>4</sup>

**FOOTNOTES**

<sup>3</sup> Allows actuator operation through the Driver Control. Includes port for the Attendant Switch.

**NOTE: All specifications and dimensions are approximate.**

## Invacare® Formula CG Tilt for TDX Spree Base Accessories and Options

CTMT-FJ Ctr Mount Adj (9.25" - 12")..... 242.00 <sup>5</sup>

### CENTER PIVOT STYLE (PAIR)

90SWING 90° Swingaway Footrests..... Standard <sup>6</sup>

### FOOTBOARD

90FB Flip up Removable Footboard..... 165.00

### FOOTPLATE OPTIONS (PAIR)

1651 Composite Flip-Up Footplates ..... Standard

AT5543 Adj Angle Flip-Up Footplates..... 205.00

CTMTFB Footboard for Center Mount..... 165.00 <sup>7</sup>

### FOOTREST/LEGREST ACCESSORIES

1337 Calf Strap ..... 39.00

1600BK Heel Loops w/Ankle Straps ..... 72.00

IMPACTG Impact Guards ..... 72.00

ALPT3 3" Extension for 90° Footrests ..... 99.00

### MISCELLANEOUS ACCESSORIES

2ATP Perspex Tray – Swing away..... 545.00

### FOOTNOTES

<sup>4</sup> Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Available on low seat-to-floor only.

<sup>5</sup> Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Available on low seat-to-floor only.

<sup>6</sup> Standard with adjustable angle footplates.

<sup>7</sup> One piece footboard for Center Mount (CTMT-F) Riggings only.