

Action 3NG

GREATER DRIVING PERFORMANCE



The Action 3NG manual wheelchair offers a sleek and stylish chair for everyday use. Delivering outstanding drive performance with new lightweight options for individuals with a more active lifestyle. Highly functional yet elegantly designed the Action 3NG wheelchair is a lightweight chair offering a vast range of options.



**Making Life's
Experiences Possible**

Yes, you can.®



Action 3NG

Adaptable and modular, the Action®3NG accommodates progressive changes in user needs. It is also simple to refurbish, with a wide range of options.

With a low transfer weight the Action®3NG is easy to fold and to dismantle, making it easily transportable in a small car.

Features



Brake and Brake extension. Fold the brakehandle for easier side transfers.



Action 3NG Lever Drive.



20 rear wheel positions.



Magnesium castor attachment.

Tech Specs



Seat width
40.5cm, 45cm,
50.5cm



Seat depth
40.5cm, 45cm



Seat height
48.5cm



Backrest height
43cm



Total width
seat width + 16cm
(transit)
seat width + 19cm
(self-propel)



Total length
108cm



Max. user weight
125kg



Transport weight
10.3kg (transit)
7.1kg (self-propel)



Total weight unloaded
13.5kg (transit)
13kg (self-propel)



Turning radius
87cm



Max. safe slope
9°



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®