

# Essential Touch

**Double-sided castellated foam mattress**

A practical support surface offering  
comfort, and pressure redistribution.



**A cost-effective solution**

For pressure injury prevention



**Yes, you can.®**

# Affordable quality mattress



## Pressure Redistribution

- › Keyhole cut cells move independently to distribute weight evenly.
- › Castellated foam with independent cells offer greater immersion and envelopment.
- › Larger outer castellations and smaller inner castellations where pressure is more prominent.



## Shear & Friction Reduction

- › Castellations help reduce shear when repositioning on the bed.
- › Multi-stretch, water-resistant, polyurethane outer cover further helps to reduce shear and friction.
- › Suitable for multiple care environments.



## Comfort

- › Castellations provide greater comfort and immersion.
- › Double-sided castellated foam allows the mattress to be turned or rotated inside the cover for increased product longevity.
- › Cover with two reinforced handles and zipped around three sides for easy handling and removal for laundering.

## Technical data ▼

For more information about this product, including the user manual, please visit our website: [www.invacare.com](http://www.invacare.com)



Sizes available ▼

2000 × 880 × 150 mm  
2000 × 1030 × 150 mm



Max. User Weight ▼

150 kg



Max. Cleaning Temperature ▼

95°

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177 (ignition Source 5)

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

### FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

© 2025 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Essential Touch - APAC - 07/2025



**Yes, you can.®**