

Invacare® Matrx® | Hardware Options



When clearance or space is an issue.

Alternative hardware enables backs **PB, PB Deep, Elite** and **Elite Deep** to be mounted when, clearance or space is an issue or for heavy duty use.

The "Easy Set" hardware comes standard with **Elite, Elite Deep, PB** and **PB Deep** backs. The "Easy Set" Heavy Duty is standard with the **Heavy Duty PB** and **PB Deep** backs.

Standard Hardware Options



1/8" mounting hardware for posture back with canes
Code: PBMH-Q



Split mounting hardware
Code: SMH



Extended mounting pins
Code: EMP



Heavy duty "Easy Set" mounting hardware
Code: HDMH



Fixed space saver mounting hardware
Code: FSSH

ADJUSTABILITY



Depth
Adjusts 2" (5cm) behind back canes



Depth
Adjusts 2" (5cm) in front of back canes



Angle
20° anterior range of back angle adjustment in front of back canes



Angle
20° posterior range of back angle adjustment behind back canes



Rotation
Back may be oriented to accommodate pelvic/trunk rotation

Hardware Options



HEAVY-DUTY EASY SET HARDWARE (PAIR)
Fits 3/4", 7/8" and 1" back canes.
Code: HDMH

EASY SET HARDWARE (PAIR)
Fits 1 1/8" back canes.
Code: PBMH-Q

EASY SET HARDWARE (PAIR)
Fits 1 1/4" back canes.
Code: PBMX-X



FIXED SPACE SAVER MOUNTING HARDWARE (PAIR)
Permanent mount recommended for Power Wheelchair application or when portability is not required.

For Matrx PB & Elite Standard or Heavy-Duty Backs
Code: FSSH

For Formula Tilt
Code: FSSH-FT

For Formula Tilt/Recline
Code: FSSH-FTR



ESR MOUNTING HARDWARE
Allows heavy-duty PB & Elite Backs to be mounted to Motion Concepts TRX Power Positioning with ESR.
Code: ESRMH-HD



REMOVABLE MOUNTING KIT FOR POWER TILT/RECLINE
Available for Matrx PB & Elite Heavy Duty Backs exclusive of PB Deep Heavy Duty.

For Permobil Tilt/Recline
Code: PERMMH

For Quantum Tilt/Recline
Code: QUANMH


For Quickie Tilt/Recline
Code: QUIKMH



MX1 Hardware Options



MONO MOUNTING HARDWARE
Includes hardware for 3/4" cane diameter
Codes: RGDMT-9
RGDMT-12
RGDMT-16




DUAL MOUNT MOUNTING HARDWARE (PAIR)
Includes hardware for 1" cane diameter
Code: FBCMT



QUICK RELEASE DUAL MOUNTING HARDWARE (PAIR)
Includes hardware for 7.8", 3/4" and 1" cane diameter
Code: RBCMT



BAND CLAMP ONLY (EACH)
Clamp for 7/8" cane diameter
Code: BCO-78
Clamp for 3/4" cane diameter
Code: BCO-34
Clamp for 1" cane diameter
Code: BCO-1



SUPERLITE MOUNTING HARDWARE (PAIR)
Includes hardware for 1" cane diameter
Code: SLMT-01
SUPERLITE MOUNTING HARDWARE (PAIR) FOR LASHER WHEELCHAIR
Includes hardware for 1 1/4" cane diameter
Code: SLMT-LW
Clamp for 1" cane diameter
Code: BCO-1

Accessories



FIXED
2.5" (6.4 cms) of vertical adjustment on back pan. Mount inside back shell without interfering with depth/angle adjustment. Available in standard 5" h x 7.5" l and paediatric 3" h x 5.5" l sizes.
Code: PBL-S-L or PBL-S-R



OFFSET FIXED
Provides an additional 3" (7.6 cms) of vertical position (higher or lower) relative to mounting slots
Available 5" h x 7.5" l.
Code: PBOL-L or PBOL-R



LATERAL PADS
Lateral Pad sizes:
(Not available on Deep models).
Code:
XS: 3.25" high x 4.25" long
S: 3.5" high x 5.25" long
M: 3.75" high x 5.75" long
L: 4.25" high x 6.5" long



SWING-AWAY
2.5" (6.4 cms) of vertical adjustment on back pan
Mount inside back shell without interfering w/depth/angle adjustment
Easy push button release 4 pad sizes.
Code: SALEB-SS-X
SS= Size XS,S,M,L
X= L Left R Right



OFFSET SWING-AWAY
Provides an additional 3" (7.6 cms) of vertical position (higher or lower). Mount inside back shell without interfering w/depth/angle adjustment. Easy push button release, 4 pad sizes.
Code: OSALEB-SS-X
SS= Size XS,S,M,L
X= L Left R Right

Accessories Continued



WIPEABLE SLIPCOVER
Moisture resistant, wipeable cover with a zipper. Available from Matrix PB or Elite Backs
Code: WSC
Ordered with Back
Code: WSCU



ELITE BACK COVER - SUPERSOFT FOAM
1" soft foam over 1" firm HR foam
Code: EBSS



EXTENDED MOUNTING PINS
Allows Matrix PB or Elite Backs to be mounted on wheelchairs that are 2" wider than the back model width. Only for backs 20" wide and 300lbs weight capacity
Code: EMP



ELITE BACK CHEST HARNESS INTERFACE PLATES
Chest harness not included - mounting kit only
Code: CHIEB



POSITIONING HARNESS FOR PB
Chest harness not included - mounting kit only
Code: CHIPB



MATRIX MINI CHEST HARNESS INTERFACE KIT
Chest harness not included - mounting kit only
Code: CHIMI



SMALL HEADREST ADAPTER FOR ELITE BACK
For Attachment of headrests to Elite 12" wide backs
Code: EHRA-S

MEDIUM HEADREST ADAPTER FOR ELITE BACK
For Attachment of headrests to Elite 13" - 17" wide backs
Code: EHRA-M

LARGE HEADREST ADAPTOR ELITE BACK
For Attachment of headrests to Elite 18" - 22" wide backs
Code: EHRA-L



PRIVACY FLAPS
X-Small: 10"W x 16"H - fits 12"-14" wide backs
Code: PFLAP-XS
Small: 13"W x 16"H - fits 15"-16" wide backs
Code: PFLAP-S
Medium: 15"W x 16"H - fits 17"-18" wide backs
Code: PFLAP-M
Large: 17"W x 16"H - fits 19"-20" wide backs
Code: PFLAP-L
X-Large: 19"W x 16"H - fits 21"-22" wide backs
Code: PFLAP-XL
XX-Large: 21"W x 16"H - fits 23"-24" wide backs
Code: PFLAP-XXL