

MX2 Back

PEAK PERFORMANCE AND LOOKS GREAT



The Invacare® Matrix® MX2 Back is the latest in back support technology. Designed for active users who demand the best, it combines aesthetic appeal with superior support and function. Lightweight carbon fibre is specially engineered to provide unmatched stability for peak performance.



MX2 Back



Cover:

- Breathable Spacer Mesh fabric, 9" and higher backs include leatherette side panels
- Cover and foam overlap back shell fo privacy, protection and comfort
- Includes removable lumbar positioning pad
- Optional privacy flap available

Carbon Fiber Back Shell:

- Exceptional durability and aesthetic appeal
- Very light weight
- 3" Contour to provide stability to hips/trunk without restricting upper body mobility

Hardware:

- Choice of Fixed Mount or Quick Release Hardware. Hardware does not come standard with the MX2 back, this will need to be selected.

Accessories

Fixed Lateral pads



Carbon Fibre Hardware



Fixed Mount

Quick Release




Privacy Flap



Tech Specs

MX2 Back



Height
6, 9, 12, 16"



Width
14 - 18"



Depth adjustment
2" Fixed HW
4" Quick-Release HW



Angle adjustment
30°



Product weight
from 680g (1.5lb)



Weight Capacity
114kg (250lb)



Limited Warranty
2 years



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®