

PB Posture Back

EXCEPTIONAL SUPPORT AND COMFORT



The Invacare® Matrix® PB Posture Backs are contoured to provide midline support and positioning. The extensive range of sizes provides exceptional support and comfort for a variety of customers.





PB Posture Back

Foam:

- Molded foam wraps around shell for maximum comfort and protection
- Includes removable positioning sacral/lumbar support pad

Outer cover:

- Moisture resistant and breathable
- Optional privacy flap is available

Back shell:

- Contoured shell provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes 5" contour depth and built-in plate for headrest attachment

Hardware:

- "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair

Accessories

Chest Harness Bracket



Privacy Flap



Split



Heavy-Duty Easy Set



Extended Range





Fixed Space Saver



Mounting Pins



Tech Specs

	 Height 12", 16", 20"	 Width 14", 16", 17", 18", 20"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40° range**	 Product weight From 1.2kg (2.8lb)	 Weight Capacity 136kg (300lb)	 Limited Warranty 2 years
PB	12", 16", 20"	14", 16", 17", 18", 20"	2" forwards, 2" backwards	40° range**	From 1.2kg (2.8lb)	136kg (300lb)	2 years
PB Heavy Duty	16", 20"	21", 23", 25"-30"	4"*** 2" forwards, 2" backwards	40° range**	from 2.7kg (6lb)	227kg (500lb)	2 years

** Limitation of depth/angle adjustment if used on chairs narrower than model width



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®