

Elite TR Back

IMPROVED IMMERSION AND COMFORT



The Invacare® Matrix® Elite TR is specially designed for use in manual or power tilt only or tilt/recline applications. It includes an extra layer of SuperSoft foam for improved immersion and comfort, and wider support across the shoulder area.



matrix®
SEATING SERIES

Elite TR Back



Outer cover:

- Oversized cover/foam featuring SuperSoft HR foam layered over HR foam for comfort and protection.
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

Back shell:

- Stylish, lightweight and breathable aluminium back shell with 3" contour depth
- Wider support in shoulder area, and slotted for lateral support hardware attachment.

Hardware:

- "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair

Accessories

Chest Harness Bracket



Privacy Flap



Split



Heavy-Duty Easy Set



Extended Range



Fixed Space Saver



Mounting Pins

















Lateral Pads



Swing Away Lateral Pads



Tech Specs

Elite TR	 Height 16", 18", 20"	 Width 16" - 20"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40°	 Product weight from 1.6kg (3.5lb)	 Weight Capacity 136kg (300lb)	 Limited Warranty 2 years
Elite TR Heavy Duty	 Height 20"	 Width 21" - 24"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40°	 Product weight from 3.9kg (8.5lb)	 Weight Capacity 227kg (500lb)	 Limited Warranty 2 years



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®