Action 3 Jnr

CREATIVE THINKING

Invacare's Action³ Junior is a lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years. Action³ Junior has been developed to match the individual needs of the child and can grow as they grow. The Action³ Junior comes in attractive, modern styles, designed to appeal to kids.





Action 3 Jnr

As children grow and develop, a key feature of a children's wheelchair is that it can be adjusted to accommodate the child's individual growth and needs. With the Invacare Action³ Junior, the width of the seat and backrest can be adjusted, as can the height of the arm and footrests.

A wide range of options are available on the Action³ Junior to meet every child's needs in terms of comfort and positioning.

Features and Options



Seat-width adjustment

The seat width of the chassis can be adjusted with an innovative telescopic crossbar over 4 seat positions from 205 to 380 mm.



Upholstery adjustment

As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps.



Front stabilizer kit

A front stabilizer kit is available to increase forward stability for active children who are likely to climb all over the wheelchair and be especially energetic when propelling the wheelchair.



Adjustable seat height and balance

For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions

Tech Specs



Seat width



Seat depth



Seat height 43.5cm (self-propel) 46cm (transit)



Backrest height 37, 40cm



Total width seat width + 17cm



96cm



Max. user weight 80kg



Transport weight 8.5kg (self-propel) 12kg (transit)



unloaded 15kg



Turning radius 75cm



Max. safe slope



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com Thailand: Phone +66 0 2108 8555 or visit www.invacare.com

